


# Journey™ Series – AntiThrust

Pressure Relieving



If an individual has a fixed kyphotic deformity, the wheelchair back should be reclined to increase the individual's visual field.

# Journey<sup>TM</sup> Series – AntiThrust

Wheelchair Cushion

## Product Specifications

## HCPC E2601 (General Use) 3" Only

Wound Risk Level:	Medium
Wound Stage:	Up to Stage II
Cushion Thickness / Weight Limit:	3" / 275lbs      4" / 300lbs

## Product Description

The Journey AntiThrust cushion is designed to prevent ischials from thrusting forward. The Ischial well is filled with Tru-Relief™ visco elastic memory foam which comfortably reduces movement in mild sliders. The unique flat-top design facilitates foot propulsion of the wheelchair while minimizing the risk for hamstring pain or tightness.

## Models and Sizes Available ( WIDTH x DEPTH x HEIGHT )

### Journey AntiThrust 3" (HCPC E2601)

Model #	Size
ST3 – 16 x 16	16" x 16" x 3"
ST3 – 16 x 18	16" x 18" x 3"
ST3 – 18 x 16	18" x 16" x 3"
ST3 – 18 x 18	18" x 18" x 3"
ST3 – 20 x 18	20" x 18" x 3"

### Journey AntiThrust 4"

Model #	Size
ST4 – 16 x 16	16" x 16" x 4"
ST4 – 16 x 18	16" x 18" x 4"
ST4 – 18 x 16	18" x 16" x 4"
ST4 – 18 x 18	18" x 18" x 4"
ST4 – 20 x 18	20" x 18" x 4"

Comfort • Safety • Mobility

Toll-Free 1(866)330-5336

**Keen**<sup>®</sup>

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

**Seating 63**