



Journey™ Series – Hard Back Insert

Hard Back Cushion Insert



The correct back height and back style will help to maintain proper posture and trunk support.

To check:

An individual should be able to insert 4 to 5 fingers between the axilla and the top of the back upholstery.

Journey™ Series – Hard Back Insert

Hard Back Cushion Insert

Product Specifications

Positioning Needs: Medium

Thickness: 1"

Additional Lumbar Depth: 1"

Product Description

The Journey Hard Back Insert is designed to provide posterior lateral pelvic stabilization along with lumbar support. With no hardware components, this back support is a quick solution to manual wheelchairs with a slinged back.

The Journey Hard Back Insert attaches to the wheelchair with four adjustable anchor straps. The aluminum core support is encased in two layers of foam.

Models and Sizes Available

Model #	Size (WIDTH x HEIGHT)
WBH – 18 x 18	18" x 18"

Comfort • Safety • Mobility

Toll-Free 1(866)330-5336

Keen

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

Seating 41