


# Journey™ Series – Wedge

Wheelchair Cushion



If the individual uses their feet for wheelchair propulsion, be sure to not use too low of a seat-to-floor height. The individual will have more problems when the wheelchair or pelvis is too low.

For example:

- It could potentially aggravate hamstring tightness and pain.
- It can lead to knee and/or hip contractures.

## Product Specifications

**HCPC E2601 (General Use)**

Wound Risk Level:	Low
Wound Stage:	Up to Stage I
Cushion Thickness / Weight Limit:	4" Front – 2" Back / 300lbs

## Product Description

The Journey Wedge cushion is designed to offer comfort and support while reducing sacral sliding. Tru-Relief<sup>TM</sup> visco elastic memory foam conforms to bony prominences while wedge design keeps an individual's pelvis in neutral pelvic alignment. The Journey Wedge reduces extensor tone by increasing hip and knee flexion.

## Models and Sizes Available ( WIDTH x DEPTH x HEIGHT )

Model #	Size
SW4 – 16 x 16	16" x 16" x 4" front to 2" back
SW4 – 16 x 18	16" x 18" x 4" front to 2" back
SW4 – 18 x 16	18" x 16" x 4" front to 2" back
SW4 – 18 x 18	18" x 18" x 4" front to 2" back
SW4 – 20 x 18	20" x 18" x 4" front to 2" back

**Comfort • Safety • Mobility**

**Toll-Free 1(866)330-5336**



Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

**Seating 65**