


Liberty™ Series – Bed Wedges

Body Wedge Cushion



“Balancing against resistance promotes stability. The side lying position may be used to initiate rolling toward supine or prone where the individual cannot initiate against the resistance.”

(Voss, Ionta, Myers 1985. Proprioceptive Neuromuscular Facilitation. Philadelphia: Harper + Row, Publishers. p.227)

Liberty™ Series – Bed Wedges

Body Wedge Cushion

Product Specifications

Purpose: Positioning

Wedge Size: 12"W x 36"L x 7"H

Product Description

The Liberty Bed Wedges are designed to offer support, positioning and stabilization of the trunk and extremities. Use in various positions as needed. The top layer is comprised of visco elastic foam, while the base is comprised of high resiliency foam.

Models and Sizes Available

Model #	Size (WIDTH x LENGTH x HEIGHT)
LBW – M	12" x 36" x 7"

Comfort • Safety • Mobility

Toll-Free 1(866)330-5336

Keen

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

Bed Safety 23