




# Liberty™ Series – Leg Arrow

Leg Flotation Cushion



4-Way Stretch Wicking Layer



Waterproof  
Breathable  
Layer



Abduction of the lower extremities can help reduce abnormal tone, scissoring of the legs, excessive internal rotation, and muscle spasms. Sleeping with a pillow or cushion between the legs can help minimize lower back pain.

# Liberty<sup>TM</sup> Series – Leg Arrow

Leg Flotation Cushion

## Product Specifications

Individual Risk Level: Prevention

## Product Description

The Liberty Leg Arrow cushion is designed to offer comfort, positioning and support of the lower extremities. Tru-Relief<sup>TM</sup> visco elastic memory foam conforms to the individual, allowing abduction and pressure redistribution from the knee to the ankle. The unique design allows the cushion to be used with individuals lying on their side, windswept in the bed or wheelchair, isolating one heel for flotation or as an abductor for individuals lying on their back.

Ideal for individuals recovering from hip surgeries, at high risk for pressure wounds, or individuals in need of corrected leg positioning. The Liberty Leg Arrow can be used with all Liberty Heel Products for total lower extremity pressure redistribution.

## Models and Sizes Available

Model #	Size
LLA – 24	One Size Fits All

**Comfort • Safety • Mobility**

**Toll-Free 1(866)330-5336**

**Keen**<sup>®</sup>

Copyright © 2007, Keen Mobility Company. All rights reserved. Rev 070312. All trademarks are property of their respective owner.

**Bed Safety 35**